## VACUSED™ for LINEAR

## PRESENTACION / CONTENTS / PRESENTATION

REF

VACUSED

5 x 100

Sólo para uso diagnóstico *in vitro*  
 Only for *in vitro* diagnostic use  
 Uniquement pour usage diagnostique *in vitro*

## VACUSED

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## INTENDED USE

VACUSED tubes for the collection of venous blood used in Erythrocyte Sedimentation Rate (ESR) procedures on Linear reading systems.

## PRINCIPLE

In the process of blood collection, these tubes work on the principle of negative pressure, collecting the blood by linking the needle to the vein.

## COMPOSITION

- 5 Pouches (Ref 5230010 1x100) Sterile plastic vacuum tubes
- Vacuum draw: 1.5 mL
- Anticoagulant: 4 NC- Citrate 3 Na

## STORAGE AND STABILITY

Store at 5-25°C.

Store the tubes in the sealed pouches until their use.

Stable up to the expiration date stated on the label. Once you take out the necessary tubes, ensure proper closure to keep tubes humidity and therefore their optimum conditions.

**Warning:** The tubes out of the pouches have a stability of 4 months (PAO).

**Do not freeze. Do not expose to excessive heat.** Avoid prolonged exposure of opened tubes to light.

*The vacuum stability depends on tube intrinsic factors as well as on tube's correct manipulation.*

## MATERIALS REQUIRED (Not provided)

- Needle, tourniquet and antiseptic
- ERIX Ref: 531005 Orbital mixer

- Visual determination: VISUAL-SED rack. Ref: 5310010
- Automatic determination: IRIA, THERMA, LENA y ERILINE

## PROCEDURE

Leave the tubes outside the pouch for 10 minutes at room temperature before use. The sample of whole blood with anticoagulant should be mixed by slowly inverting the tube during approximately 5 minutes before starting the assay. Follow the instructions of the analyzer.

## CAUSES OF ERROR

- Dirty material (re-use of disposable material)
- Haemolysed samples
- Bubbles in samples
- Insufficient sample mixing

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