

Ferritin - Turbidimetric

<p>REF 3140005 1 x 50 mL CONTENTS R1. Reagent 1 x 40 mL R2. Reagent 1 x 10 mL CAL. 1 x 3 mL</p>	<p>REF 3140010 2 x 50 mL CONTENTS R1. Reagent 2 x 40 mL R2. Reagent 2 x 10 mL CAL. 1 x 3 mL</p>
For <i>in vitro</i> diagnostic use only	

Ferritin - Turbidimetric

Latex Turbidimetry

PRINCIPLE

The latex particles coated with anti human ferritin are agglutinated when they react with samples that contain ferritin. The latex particles agglutination is proportional to the concentration of the ferritin in the sample and can be measured by turbidimetry.^{1,2}

REAGENTS COMPOSITION

- R1** **Diluent.** Glycine buffer, 20 mmol/L, pH 8.5.
- R2** **Latex.** Latex particles coated with polyclonal anti-human ferritin antibodies, pH 8.2.
- CAL** **Calibrator.** Ref 3940005 1x3 mL. Human ferritin. Ferritin concentration is stated on the label vial and it is trazable to the 3rd International Reference Material for Ferritin, 94/572 WHO (NIBSC).

Precautions: The reagents contain sodium azide 0.95 g/L. Avoid any contact with skin or mucous. The reagents from human donors have been given negative results to anti-HIV 1/2, HBsAg and anti-HCV. Handle cautiously is recommended.

REAGENT PREPARATION

- R1** Ready to use.
- R2** Ready to use. Shake gently the vial before use.
- CAL** Ready to use.

Calibration curve: Prepare dilutions of the Calibrator using NaCl 9 g/L as diluent. Multiply the concentration of the Calibrator by the corresponding factor indicated in the table below to obtain the ferritin concentration of each point of the curve.

Dilution	1	2	3	4	5
Ferritin-CAL (µL)	--	25	50	75	100
NaCl 9 g/L (µL)	100	75	50	25	--
Factor	0.0	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

- The reagents will remain stable until the expiration date printed on the label, when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use the reagents after the expiration date.
- Reagent deterioration: Presence of particles, turbidity and increment of blank reagent.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing. Hemolyzed or contaminated samples are not suitable for testing.

MATERIAL REQUIRED

- Thermostatic bath at 37°C.
- Spectrophotometer or photometer thermostatable at 37°C with a 650 ± 20 nm filter.

PROCEDURE

Preliminary Procedure

Prewarm the reagents and the photometer (cuvette holder) to 37°C.

Analytic Procedure

- Using distilled water zero the instrument at 650 nm.
- Pipette into a cuvette:

Diluent: R1	0.8 mL
Sample/ Calibrator/ Water (Blank)	100 µL
Latex: R2	0.2 mL

- Mix well and record the absorbance immediatelly (A₁) and after 8 minutes (A₂) of the reagent R2 addition.



Calculation

Calculate the absorbance difference (A_2-A_1) of each point of the calibration curve and plot the values obtained against the ferritin concentration of each calibrator dilution. Ferritin concentration in the sample is calculated by interpolation of its (A_2-A_1) in the calibration curve.

QUALITY CONTROLS

Control sera are recommended to monitor the performance of manual and automated assay procedures. It is recommended to use Plasma Protein Control N-I (ref: 3915010) and N-II (ref: 3915015). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES^{3,7}

Children: 7 – 140 µg/L

Men: 20 – 250 µg/L

Women: 20 – 200 µg/L

Each laboratory should establish its own reference range.

CLINICAL SIGNIFICANCE³⁻⁷

Ferritin is the major iron storage compound in the body and is considered one of the most reliable indicators of iron status of patients.

A clinical evaluation of serum ferritin is an index of iron stores.

Whereas low serum concentrations of ferritin are always indicative of an iron deficiency, elevated concentrations can occur for variety of reasons. Thus, although elevated concentrations often indicate an excessive iron intake, they are also caused by liver disease, chronic inflammation and malignancies. Pregnant women, blood donors, hemodialysis patients, adolescents and children are groups particularly at risk. Plasma ferritin is also increased in patients with hemosiderosis or hemochromatosis.

ANALYTICAL PERFORMANCE

- **Linearity limit:** Up to 300 µg/L, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in CINA 9 g/L and retested again.
- **Detection limit:** Values less than 3.0 µg/L give non-reproducible results.
- **Analytical sensitivity:** 2.07 mA / µg/L
- **Prozone effect:** Up to 4000 µg/L

Precision:

	Mean (µg/L)	CV (%)
Intra-assay N = 10	65 178	3.56 1.87
Inter-assay N = 10	65 178	5.16 2.9

- **Accuracy:** Results obtained with this reagents did not show systematic differences when compared with commercial reagents of similar characteristics. Studies of comparison are available on request.
- **Interferences:** Bilirubin (20 mg/dL), hemoglobin (10 g/L) and rheumatoid factors (600 IU/mL), do not interfere. Lipemia interferes. Other substances may interfere⁸.

NOTES

1. Calibrator dilutions in plastic tubes should be avoided. Ferritin antigen may coat to the walls of plastic tubes.
2. This method may be used with different instruments. Any application to an instrument should be validated to demonstrate that results meets the performance characteristics of the method. It is recommended to validate periodically the instrument. Contact to the distributor for any question on the application method.
3. The linearity limit depends on the sample/reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
4. Heterophilic antibodies in human serum can react with reagent antibodies, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed. Additional information may be required for diagnosis.
5. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Newman DJ, Henneberry H, Price CP. *Ann Clin Biochem* 29: 122-42 (1992).
2. Bernard A, Lauwerys R. *J Immunol Methods*. 71:141-147 (1984).
3. Wiedermann, Jonetz-Mentzel. *Eur J Clin Chem Clin Biochem* 31: 453 - 457 (1993).
4. Worwood M. *Blood Reviews*. 4 : 259 - 269 (1990).
5. Lipschitz D, Cook JD, Finch CA. *The New England J of Med* 290: 1213-1216 (1974).
6. Mazza J, Barr RM, McDonald JWD. *Can Med Assoc J* 119: 884-886 (1978).
7. Tietz Textbook of Clinical Chemistry, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
8. Young DS. *Effects of drugs on clinical laboratory tests*. 3th ed. AACCC Press (1997).

T3140-2/1111
R1.ing

