

LISA SHAKE

Agitator – incubator ELISA

Moduri de lucru	Agitator
	Incubator
	Agitator si Incubator
Controlul temperaturii	
Domeniul de temperatura Rezolutia	de la 30 ⁰ C la 45 ⁰ C 1 ⁰ C
Timpul de incubare	de la 1 la 999 minute
Agitator	
Frecventa Amplitudinea	de la 400 la 700 RPM 2 mm
Pozitia de exploatare	Pe o suprafata plana, rigida si lipsita de vibratii
Conditii de exploatare	
Temperatura Umiditatea relativa	de la + 18 ⁰ C la 32 ⁰ C pana la 80%, fara condensare
Conditii de depozitare	
Temperatura Umiditatea relativa	de la -10 ⁰ C la 35 ⁰ C pana la 80%, fara condensare
Carcasa	ABS Fire retardant
Dimensiuni (cm)	28 X 26 X 14 (l X b X h)
Greutate	2.5 Kg (Approx)

