

Auto-coding Element™

Blood Glucose Monitoring System




- ☑ Domeniul de testare: 10-600 mg/dL
- ☑ Dimensiuni: 76 x 56 x 18 (mm)
- ☑ Greutate: 45g (incluzand baterii)
- ☑ Viata bateriei: 5000 teste
- ☑ Port comunicare PC

Auto
Coding

3
sec

0.3µL



- ☑ SISTEM DE AUTO-CODARE A STRIPURILOR
- ☑ PROBA DE DOAR 0.3 µL
- ☑ REZULTAT IN 3 SECUNDE
- ☑ MARCAREA CONDITIEI TESTULUI (INAINTE DE MASA, DUPA MASA, DUPA EXERCITIU FIZIC)
- ☑ LOCURI ALTERNATIVE DE TESATRE (BRAT, ANTEBRAT, PALMA, COAPSA, GAMBA)
- ☑ MEMORIE PENTRU 365 TESTE
- ☑ CODURI PANA LA 100 DE UTILIZATORI

BALMED s.r.l.

Bd. Burebista Nr. 1
Bl. D15, Sc.4, Ap. 125
031106 BUCURESTI
Tel.fax:021/327.52.69

www.balmed.ro

Diabetes Management Technology



infopia Co., Ltd.