

Lampa de fototerapie M304





Diagrama spectrala a lampii



EN ISO 9001:2000



| MODEL | M304 |
|--|--|
| Sursa de lumina de fototerapie | Lampa fluorescenta |
| Numarul de lampi albastre | 4 |
| Durata de viata a lampilor albastre | 2.000 ore |
| Tipul de lampa albastra | Philips TL 20 W/ 52 – Lampa de presiune scazuta cu vapori de mercur, cu emisie la 450 nm (UV-B) |
| Numarul de lampi lumina alba | 1 – lampa florescenta |
| Cronometru | Programabil, 24 ore |
| Indicator total ore functionare | 0-99.999 ore |
| Puterea lampii lumina alba | 20W |
| Puterea lampilor albastre | 4 x 20 W = 80 W |
| Metode de lucru | 1) direct ; 2) cu cronometru |
| Puterea absorbita | 100 W |
| Siguranta fuzibila | 1 A |
| Carcasa | Otel emailat |
| Materialul suportului | Otel emailat |
| Dimensiunile unitatii de fototerapie | (L x i x l) 65 x 13.5 x 32 cm |
| Dimensiuni generale | 78 x 200 x 60 cm (la inaltimea maxima) |
| Greutatea unitatii de fototerapie | 9,9 Kg |
| Greutatea totala neta | 27,7 Kg |
| Manual de utilizare | in limba romana si engleza |
| ELEKTRO-MAG A.Ş İkitelli O.S.B Demirciler sitesi B7 Blok No:153 İkitelli/İSTANBUL Tel : 0212 549 55 25 Fax: 0212 549 52 09 e-mail: info@elektro-mag.com.tr www.elektro-mag.com.tr | in Romania prin: BALMED s.r.l. Bd. Burebista Nr. 1 Bl. D15, Sc.4, Ap. 125 031106 BUCURESTI Tel/fax. 021 327 52 69 Tel/fax. 021 327 09 21 e-mail: office@balmed.ro www.balmed.ro |