


ACCESORII ANALIZOARE VSH AUTOMATE

LINEAR

Cromatest

clonatest

Control pe doua nivele




ERIX – mixer orbital



Viteza : fixa
 Nr. de rotatii: 6.6 rot/min
 Probe ce pot fi mixate: Tuburi nevacuumate
 ori vacuumate
Incarcarea: 20 tuburi
 Timp de mixare recomandat: min. 5 min
 Temperatura de mixat: Temperatura camerei
 Nr. de rotatii: 6.6 rot/min
 Caracteristici: 2 pozitii (inclinatii)
 Cerinte ambientale:
 Temperatura: camerei
 Raza platanului: 15.5 cm
 Masa: 2.5 kg

LENA TP - imprimanta



Imprimanta pe hartie termica dedicata
 analizoarelor IRIA si LENA
 Conexiune prin port serial RS232
 Dimensiuni 11 x 14 x 7 cm
 Masa : 410 g cu rola de hartie


Rezultate IRIA

POS 1	9	30
POS 2	85	137
POS 3	78	129
POS 4	16	45
POS 5	74	126


Rezultate pacienti LENA NE, THERMA NE

```
List Ps R1h R2h ID_Code
3 9 10 .....
4 10 52 .....
```

```
List Ps R1h R2h ID_Code
9 10 123456789111315
```




```
List Ps R1h R2h ID_Code
10 52 100000000000016
```



Rezultate QC LENA NE, THERMA NE

```
QC2
55 08/05/16
54 08/05/15
54 08/05/15
56 08/05/15
57 07/05/15
54 07/05/15
55 07/05/14
57 06/05/14
57 06/05/14
56 05/05/14
```

QC2



BALMED s.r.l.

Bd. Burebista Nr. 1
 Bl. D15, Sc.4, Ap. 125,
 031106 BUCURESTI
 Tel.fax:021/327.52.69
 www.balmed.ro